# OHST

medical technology

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**Grünauer Fenn 3**

**14712 Rathenow**

**Deutschland / Germany / Allemagne /**

**Germania / Alemania / Duitsland /**

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## IMPLANT

### Biathlon Hip Stem

Before using the product, the user is under obligation to carefully study and comply with the following recommendations and information along with the information specific to the product.

The party introducing this product into circulation accepts no liability for direct or consequential damage resulting from careless use or handling, particularly noncompliance with the following user instructions or improper care or maintenance.

These implants may be used only by physicians with appropriate experience, practice and skills in hip arthroplasty. Familiarity with the surgical technique recommended for this system and its diligent application are indispensable in order to achieve the best possible result.

#### 1. Product description and implant materials

The Biathlon Hip Stem is fixed cementless, consists of ISO 5832-3 Ti6Al4V and has a proximal TPS coating. In addition, the cementless variation is available with an optional bioactive calcium phosphate coating. There are 10 sizes obtainable in each case with two different offsets (STD and LAT).

Product, packaging contents and materials used are specified on the product label. The implant must be implanted using a suitable surgical technique familiar to the surgeon. In this regard, attention must be given to the explanations concerning the particular surgical technique.

#### 1.1 Implant overview

Name	Material	Reference number
Biathlon Hip Stem STD Sz. 1 cementless, TPS	ISO 5832-3 Ti6Al4V	367-1038
Biathlon Hip Stem STD Sz. 2 cementless, TPS	ISO 5832-3 Ti6Al4V	367-1039
Biathlon Hip Stem STD Sz. 3 cementless, TPS	ISO 5832-3 Ti6Al4V	367-1040
Biathlon Hip Stem STD Sz. 4 cementless, TPS	ISO 5832-3 Ti6Al4V	367-1041
Biathlon Hip Stem STD Sz. 5 cementless, TPS	ISO 5832-3 Ti6Al4V	367-1042
Biathlon Hip Stem STD Sz. 6 cementless, TPS	ISO 5832-3 Ti6Al4V	367-1043
Biathlon Hip Stem STD Sz. 7 cementless, TPS	ISO 5832-3 Ti6Al4V	367-1044
Biathlon Hip Stem STD Sz. 8 cementless, TPS	ISO 5832-3 Ti6Al4V	367-1045
Biathlon Hip Stem STD Sz. 9 cementless, TPS	ISO 5832-3 Ti6Al4V	367-1046
Biathlon Hip Stem STD Sz. 10 cementless, TPS	ISO 5832-3 Ti6Al4V	367-1047
Biathlon Hip Stem LAT Sz. 1 cementless, TPS	ISO 5832-3 Ti6Al4V	367-1067
Biathlon Hip Stem LAT Sz. 2 cementless, TPS	ISO 5832-3 Ti6Al4V	367-1068
Biathlon Hip Stem LAT Sz. 3 cementless, TPS	ISO 5832-3 Ti6Al4V	367-1069
Biathlon Hip Stem LAT Sz. 4 cementless, TPS	ISO 5832-3 Ti6Al4V	367-1070
Biathlon Hip Stem LAT Sz. 5 cementless, TPS	ISO 5832-3 Ti6Al4V	367-1071
Biathlon Hip Stem LAT Sz. 6 cementless, TPS	ISO 5832-3 Ti6Al4V	367-1072

<b>Name</b>	<b>Material</b>	<b>Reference number</b>
Biathlon Hip Stem LAT Sz. 7 cementless, TPS	ISO 5832-3 Ti6Al4V	367-1073
Biathlon Hip Stem LAT Sz. 8 cementless, TPS	ISO 5832-3 Ti6Al4V	367-1074
Biathlon Hip Stem LAT Sz. 9 cementless, TPS	ISO 5832-3 Ti6Al4V	367-1075
Biathlon Hip Stem LAT Sz. 10 cementless, TPS	ISO 5832-3 Ti6Al4V	367-1076
Biathlon Hip Stem STD Sz. 1 cementless, TPS+CaP	ISO 5832-3 Ti6Al4V	367-1087
Biathlon Hip Stem STD Sz. 2 cementless, TPS+CaP	ISO 5832-3 Ti6Al4V	367-1088
Biathlon Hip Stem STD Sz. 3 cementless, TPS+CaP	ISO 5832-3 Ti6Al4V	367-1089
Biathlon Hip Stem STD Sz. 4 cementless, TPS+CaP	ISO 5832-3 Ti6Al4V	367-1090
Biathlon Hip Stem STD Sz. 5 cementless, TPS+CaP	ISO 5832-3 Ti6Al4V	367-1091
Biathlon Hip Stem STD Sz. 6 cementless, TPS+CaP	ISO 5832-3 Ti6Al4V	367-1092
Biathlon Hip Stem STD Sz. 7 cementless, TPS+CaP	ISO 5832-3 Ti6Al4V	367-1093
Biathlon Hip Stem STD Sz. 8 cementless, TPS+CaP	ISO 5832-3 Ti6Al4V	367-1094
Biathlon Hip Stem STD Sz. 9 cementless, TPS+CaP	ISO 5832-3 Ti6Al4V	367-1095
Biathlon Hip Stem STD Sz. 10 cementless, TPS+CaP	ISO 5832-3 Ti6Al4V	367-1096
Biathlon Hip Stem LAT Sz. 1 cementless, TPS+CaP	ISO 5832-3 Ti6Al4V	367-1097
Biathlon Hip Stem LAT Sz. 2 cementless, TPS+CaP	ISO 5832-3 Ti6Al4V	367-1098
Biathlon Hip Stem LAT Sz. 3 cementless, TPS+CaP	ISO 5832-3 Ti6Al4V	367-1099
Biathlon Hip Stem LAT Sz. 4 cementless, TPS+CaP	ISO 5832-3 Ti6Al4V	367-1100
Biathlon Hip Stem LAT Sz. 5 cementless, TPS+CaP	ISO 5832-3 Ti6Al4V	367-1101
Biathlon Hip Stem LAT Sz. 6 cementless, TPS+CaP	ISO 5832-3 Ti6Al4V	367-1102
Biathlon Hip Stem LAT Sz. 7 cementless, TPS+CaP	ISO 5832-3 Ti6Al4V	367-1103
Biathlon Hip Stem LAT Sz. 8 cementless, TPS+CaP	ISO 5832-3 Ti6Al4V	367-1104
Biathlon Hip Stem LAT Sz. 9 cementless, TPS+CaP	ISO 5832-3 Ti6Al4V	367-1105
Biathlon Hip Stem LAT Sz. 10 cementless, TPS+CaP	ISO 5832-3 Ti6Al4V	367-1106

## 1.2 Instrument overview

The instruments of OHST Medizintechnik AG listed below must be used exclusively for implantation:

<b>Name</b>	<b>Reference number</b>
Instrument set Biathlon Hip Stem Tray I	367-289
Instrument set Biathlon Hip Stem Tray II	367-290
Hip Stem Extraction Instrument Set	206-010

## 1.3 Accessories

<b>Name</b>	<b>Reference number</b>
Surgical Technique Biathlon Hip Stem	50000509
X-ray template Biathlon Hip Stem STD cementless	367-1242
X-ray template Biathlon Hip Stem LAT cementless	367-1243
Instructions for use stem reducer	50000427

Name	Reference number
Instructions for Use Extraction Instrument Set	50000428
Implant passport	50000572

## 2. Handling

### 2.1 General information

This implant is part of a system and must only be used with the appropriate original system components. Only the instruments of the system listed above must be used for implantation. Before using the instruments the attached instructions for use (50000354) must be considered.

**Caution:** Implants must always be kept in their complete, unopened protective packaging. The packaging containing the implant must not be exposed to direct sunlight. Before inserting the implant, the packaging must be examined for damage, as this could affect sterility.

When unpacking the implant, its conformity with the designation on the packaging (art. no. / serial no. / size) must be checked. Compliance with the appropriate hygiene regulations is required during removal of the implant from the packaging. Care must be taken to protect all implant surfaces against damage, since this could be decisive for possible failure. The prosthesis must therefore not come into contact with objects which could damage its surface. Before use, every implant must be visually inspected for damage. Machining an implant can not only reduce its service life, but can also lead to immediate or subsequent failure of the prosthesis under stress. The implant must therefore neither be mechanically nor otherwise processed. Implants from damaged packaging, unsterile, contaminated, damaged or carelessly handled implants or implants subjected to unauthorized machining must not be used.

**Caution:** Implants are intended for single use only! Individual loads on the functional surfaces of an implant used for one patient modify the functional surfaces in a way that excludes any reuse. Detection of load-caused markings by visual methods only is not secured. After explantation, previous damages which exclude any reuse must therefore be assumed.

### 2.2 Authorised component combinations

The following femoral heads must be used for assembly with the implant:

Nominal Ø	Name	Neck length	Reference number
Ø22 mm	Femoral Head CoCrMo ISO 5832-12, 12/14 S / M / L	-4 / 0 / +4 mm	030-2200 to 030-2202
Ø28 mm	Femoral Head Implant steel ISO 5832-9, 12/14 S / M / L / XL / XXL	-3,5 / 0 / +3,5 / +7 / +10,5 mm	020-2800 to 020-2804
	Femoral Head CoCrMo ISO 5832-12, 12/14 S / M / L / XL / XXL	-3,5 / 0 / +3,5 / +7 / +10,5 mm	030-2800 to 030-2804
	Femoral Head Biolox® forte ISO 6474-1, 12/14 S / M / L	-3,5 / 0 / +3,5 mm	367-907 to 367-909

<b>Nominal Ø</b>	<b>Name</b>	<b>Neck length</b>	<b>Reference number</b>
	Femoral Head Biolox® delta ISO 6474-2, 12/14 S / M / L	-3,5 / 0 / +3,5 mm	367-1140 to 367-1142
	ELEC® Femoral Head ISO 6474-1, 12/14 S / M / L	-3,5 / 0 / +3,5 mm	384-001 to 384-003
	ELEC®plus Femoral Head ISO 6474-2, 12/14 S / M / L	-3,5 / 0 / +3,5 mm	110230, 110240, 110250
	ELEC®plus Femoral Head ISO 6474-2, 12/14 S / M / L	-3,5 / 0 / +3,5 mm	013-001 to 013-003
Ø32 mm	Femoral Head Implant steel ISO 5832-9, 12/14 S / M / L / XL / XXL	-4 / 0 / +4 / +8 / +12 mm	020-3200 to 020-3204
	Femoral Head CoCrMo ISO 5832-12, 12/14 S / M / L / XL / XXL	-4 / 0 / +4 / +8 / +12 mm	030-3200 to 030-3204
	Femoral Head Biolox® forte ISO 6474-1, 12/14 S / M / L	-4 / 0 / +4 mm	367-920 to 367-912
	Femoral Head Biolox® delta ISO 6474-2, 12/14 S / M / L / XL	-4 / 0 / +4 / +7 mm	367-1143 to 367-1145, 367-1149
	ELEC® Femoral Head ISO 6474-1, 12/14 S / M / L	-4 / 0 / +4 mm	384-004 to 384-006
	ELEC®plus Femoral Head ISO 6474-2, 12/14 S / M / L / XL	-4 / 0 / +4 / +7 mm	110260, 110270, 110280, 110291
	ELEC®plus Femoral Head ISO 6474-2, 12/14 S / M / L / XL	-4 / 0 / +4 / +7 mm	013-004 to 013-007
Ø36 mm	Femoral Head Biolox® forte ISO 6474-1, 12/14 S / M / L	-4 / 0 / +4 mm	367-930 to 367-932
	Femoral Head Biolox® delta ISO 6474-2, 12/14 S / M / L / XL	-4 / 0 / +4 / +8 mm	367-1146 to 367-1148, 367-1150
	ELEC® Femoral Head ISO 6474-1, 12/14 S / M / L	-4 / 0 / +4 mm	384-007 to 384-009

Nominal Ø	Name	Neck length	Reference number
	ELEC@plus Femoral Head ISO 6474-2, 12/14 S / M / L / XL	-4 / 0 / +4 / +8 mm	110300, 110310, 110320, 110330
	ELEC@plus Femoral Head ISO 6474-2, 12/14 S / M / L / XL	-4 / 0 / +4 / +8 mm	013-008 to 013-011
Ø40 mm to Ø60 mm	Fracture head size S implant steel ISO 5832-9, 12/14	-4 mm	155-140 to 155-160
	Fracture head size M implant steel ISO 5832-9, 12/14	0 mm	155-040 to 155-060
	Fracture head size L implant steel ISO 5832-9, 12/14	+4 mm	155-240 to 155-260

We only guarantee compatibility of our products in combination with our own CE-marked products and with the products we have approved for combined use, for which the competent authority has issued its approval. In this regard, please note the instructions for use of the endoprosthesis manufacturers and the combination matrix approved by OHST.

Due to reasons relating the product safety and product liability, it is prohibited to use implants manufactured by OHST Medizintechnik AG in combination with components of other manufacturer that have not been approved by OHST.

### 2.3 Information for use

The application of the implant is done without cement. A cement restrictor must be used for cemented applications. Before using the cement restrictor the corresponding instructions for use must be considered. The implant is equipped with a 12/14 cone to connect it with a femoral head. The femoral stem cone and the inner cone of the femoral head must be clean and intact when assembled. Before attaching the femoral head the cone should be cleaned carefully. The appropriate femoral head should then be attached by hand and fixed on the cone with the head insertion instrument and an appropriate hammer blow. When using ceramic femoral heads the corresponding instruction for use is also to be considered. After reduction, the definitive stability, mobility and muscle tension must be checked.

**Caution:** Please be hereby explicitly advised that, in case of an intraoperative change or revision of the femoral head, only femoral heads without a ceramic cone are to be used. This is valid irrespective of the materials used in the previous cone pairing.

**Caution:** If a ceramic component is damaged or fractured, complete revision of the prosthetic components at the earliest possible date is recommended. In this case, the use of metal femoral heads is contraindicated in revision surgery, as this may lead to serious and partly life-threatening complications. In the rare event of a fracture of the ceramic component, thorough debridement with removal of all visible ceramic particles as well as careful wound irrigation is absolutely essential during surgery.

If it should be necessary intraoperatively to remove the already inserted original prosthesis, a hip stem extractor is available. Before inserting the implant, the implant bed must be irrigated sufficiently. During

implantation, ensure that all loose particles (e. g. bone splinters, friction particles from the instruments) are removed from the prepared implant bed.

The porous coating (TPS, Bonit®, CaP, HA) of the implant surfaces and the roughened surfaces must not come in contact with clothing or other fibre-shedding materials.

**Caution:** Whenever possible, avoid contact with the BONIT®-coated implant parts. These areas should be touched only with powder-free latex gloves.

**Caution:** When using high-frequency surgical instruments (e.g. cautery knife), it must be ensured that they do not come into contact with the implants or instruments. This can cause such severe damage to the implants or instruments that failure (e.g. fracture) may ensue. If an implant has been damaged, it must not remain in the patient but needs to be replaced by a new, intact implant. Damaged instruments may only continue to be used if they can still perform their intended function without compromise.

## 2.4 Surgical technique

After opening the joint capsule and dislocating the femoral head from the cup, the femoral head is resected as determined in the preoperative planning and removed in full.

A femoral neck resection guide is optionally available for marking the resection plane. Then open the medullary cavity using the cavity chisel. The chisel has to be applied to the far lateral and dorsal side in order to facilitate the subsequent driving in of the rasps in the direction of the femoral axis. The cavity chisel should be applied according to the desired antetorsion. Fracturing of the trochanter major must be avoided. The medullary cavity opening can be opened further with the surgical awl.

Once the medullary cavity has been opened, the femoral stem is rasped. The first rasp determines the orientation of the subsequent sizes. Care must be taken to ensure correct antetorsion already before the first rasp is applied. The antetorsion can be reviewed by means of the guide rod and rasp handle and is usually 10°-15°. Begin with the starter rasp once the rasp handle has been attached to it. After that, the subsequent rasps are used in ascending order until the preoperatively determined size is reached. The sizes of the rasps conform to the sizes of the implants. The correct position of the rasp inside the femur can be verified with an image converter. If during the rasping procedure the stem size does not correspond to the preoperatively determined size and if there is a difference of two or more sizes, the alignment of the axis may be incorrect or there may be an osseous obstruction present. In this case the selected stem might be too small and would therefore not provide the required stability. On request, special rasp handles are available for the various access routes.

The rasps are designed in such a way that they can also be used for trial reductions. The leg length, the range of motion and the tension of the ligaments can be checked by attaching the various trial cones and trial heads. Remove the handle from the rasp, which remains inside the femur, and attach the trial cone to the rasp. Trial cones with two offsets are offered for trial reductions of the various stem versions. The trial cones snap into place as soon as they are in their correct position. Then attach the trial head to the trial cone by hand. Trial heads with various diameters are offered in the neck lengths S to XXL for trial reductions. After reduction, the definitive stability, mobility and muscle tension should be checked.

Insert the prosthetic device into the implant bed up to the resection plane using the stem impactor. After that, thoroughly clean and dry the prosthetic cone and attach the femoral head with the previously determined neck length (S to XXL). Then reduce the stem with the femoral head into the cup and check the range of motion, leg length and tension of the ligaments. The operation is routinely completed with the layer-by-layer wound closure.

### **3. Packaging and sterility**

Depending on the sterilisation method used, implants are packaged in a triple transparent pouch made of plastic laminated film (sterilisation by irradiation at least 25 kGy) or in a double transparent pouch made of Tyvek® (ethylene oxide sterilisation) with a carton. The instruments are supplied unsterile in protective packaging. They must be cleaned and sterilised prior to use in accordance with the respective instructions for use (50000354). The stated expiry date presumes that the packaging is intact and unopened and that the product is stored under suitable conditions.

**Caution:** The implants may not be resterilized! The reconditioning of components that have not been implanted but the packaging of which has been opened is permitted only at the manufacturer, because the components must pass through individual validated processes once again.

The outer pouch of the triple transparent pouch packaging is to be removed by the non-sterile personnel together with the carton. For the double transparent pouch packaging, only the carton is to be removed by the non-sterile personnel. The second pouch must be opened such that the sterility of the inner pouch is not compromised. The inner pouch is removed and opened by the sterile personnel. The implant must then be presented to the surgeon, who can then directly remove the sterile implant.

### **4. Preoperative planning and postoperative care**

Preoperative planning by reference to X-rays, CT data and similar data is indispensable and provides important information about suitable implants, placement, possible component combinations and enables the size of the implant to be used to be preselected. The operation may be performed only if the patient's material compatibility has been determined. Use the X-ray templates for planning the operation. These are available for all sizes in a magnification of 1.15:1. In addition, X-ray templates with a 1:1 ratio are available in digital form. Trial prostheses for checking the correct seating (where applicable) and additional implants should be available in case another size is required or the intended implant cannot be used. Recognized procedures must be used for postoperative care.

### **5. Indications**

- Advanced wear of the hip joint based on degenerative, post-traumatic or rheumatoid arthritis (coxarthrosis)
- Fracture or avascular necrosis of the femoral head
- Sequelae of earlier surgical procedures, e.g. osteosynthesis, articular reconstruction, arthrodesis, hemiarthroplasty or total hip prosthesis

### **6. Contraindications**

- Acute or chronic infections, locally or systemic
- Severe muscle, nerve or vascular diseases endangering the extremity concerned
- Missing bone substance or poor bone quality that threatens the stable fit of the prosthesis
- Hypersensitivity to the materials used
- Any underlying condition that can endanger the function of the implant
- Revision with extensive bony defects

If femoral heads with a neck (XL and XXL) are used, the range of motion is reduced by about 30° and achieves values between 80° and 100° in flexion and extension.

## **7. Risks and conditions that can impair the success of the operation**

Potential risks associated with the operation are:

- Disorders of bone metabolism (osteoporosis, osteomalacia)
- Occurrence of fissures, in rare cases fracture of the femoral bone
- Circulatory disorders of the affected limb
- Neurological disorders of the affected limb
- Muscle malfunction in the affected limb
- Overweight
- Alcohol or substance abuse
- Patient groups with mental disorders or addictions
- Pregnancy
- Growth in children and adolescents
- Anticipated extreme loading e.g. due to work and sport
- Epilepsy or other reasons for repeated trauma with an increased risk of fracture
- Joint deformities that make fixation of the implant difficult
- Weakening of the bearing structures by tumour
- High-dose ingestion of cortisone or cytostatics
- Previous or threatening infectious diseases with possible joint involvement
- Deep vein thrombosis and/or history of pulmonary embolism
- All general surgical risks

## **8. Undesirable effects**

The adverse events listed below are among the most typical and commonly occurring consequences of total hip arthroplasty:

- Change in position and loosening of the prosthesis
- Dislocation of the prosthesis
- Implant breakage
- Infection
- Venous thrombosis and pulmonary embolism
- Cardiovascular disorders
- Haematomas
- Paresthesia
- Numbness
- Swelling
- Nerve damage
- Muscle spasm
- Stiffness,
- Implant noises
- Reduced quality of life (pain, sleep disorders, ROM limitations; in particular also when lying down)
- Inflammation
- Oedem/fluids
- Metallosis

- Elevated metal ions in blood
- Coxa vara
- Osteolysis
- Heterotopic ossification
- Pseudotumours

## **9. Patient information, documentation**

The serial numbers of the implants used must be recorded in the patient's records. Appropriate labels are included with the packaging of the sterile implants.

The patient must be informed of the advantages and risks of the procedure. If the implant is regarded as the best solution for the patient, even though the contraindications described above partially apply to the patient, the patients should be informed about the possible effects of these circumstances and the anticipated risks. It should be explained to patients receiving a hip replacement that the life of the implant will depend on their weight and degree of activity. It is also recommended that the patient is informed of activities with which he can reduce the effects of these aggravating circumstances.

All the information given to the patient must be documented in writing by the surgeon. After surgery, the patient must be given an implant pass containing all necessary information concerning the implant. Adverse effects that are harmful to patients can arise during MRI investigations. Artefacts, heating of implant, induction of electrical currents and implant loosening are among the possible effects. The equipment manufacturer's instructions should be carefully studied before use. In case of doubt, comparable implants should be checked for their specific MRI suitability as part of an individual risk assessment. Patients should be informed of the risks.

## **10. Key to label symbols**

An explanation of the symbols used by OHST Medizintechnik AG can be found in the annex (p. 179).

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