OHST 
medical technology

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IMPLANT

ZEN Knee System

Before using the product, the user is under an obligation to carefully study and follow the recommendations and information below, together with the information specific to the product.

The party introducing this product into circulation accepts no liability for direct or consequential damage or injury resulting from careless use or handling, and in particular non-compliance with the following instructions for use or improper care or maintenance.

These implants may only be used by physicians with detailed knowledge, experience and the skills required for knee arthroplasty. Familiarity with the surgical technique recommended for this system and its careful application is essential in order to achieve the best possible result.

Care must be taken to always use the most up-to-date version of the instructions for use provided with the system components. The date of revision is shown on the cover sheet in the footer.

1. Product Description and Implant Materials

Knee prostheses are implants that remain in the body permanently and either partially or fully replace the damaged knee. The knee joint in its entirety is replaced by three components: the femoral component, the tibial component and the patellar component.

Irrespective of the shape and fixation manner of the individual joint components, they always perform the same functions (ISO 7207-1, ISO 21536):

Femoral component: the component of a total knee joint replacement intended to be secured to the femur to replace its articulating surfaces. These implants can be manufactured as one component or a set of components to be assembled by the user.

Tibial component: the component of a total knee joint replacement intended to be secured to the tibia to replace its articulating surfaces. These implants can be manufactured as one component or a set of components to be assembled by the user.

Patellar component: the component of a total or partial knee joint replacement which is used to replace the articulating surface of the patella. These implants can be manufactured as one component or a set of components to be assembled by the user.

The product, the package contents and the materials used are specified on the product label. The implant is to be inserted using a suitable surgical technique with which the surgeon is familiar. To this end, the explanations in the respective surgical instructions must be followed.

1.1 Overview of Implants

1.1.1 Femoral components

Femur STD:

Size	Version	Reference number		
		cemented		cementless
		ISO 5832-4 CoCrMo	ISO 5832-4 CoCrMo TiNbN	ISO 5832-4 CoCrMo TPS, BONIT®
1	right	482-002	482-122	482-014
2	right	482-004	482-124	482-016
3	right	482-006	482-126	482-018
4	right	482-008	482-128	482-020
5	right	482-010	482-130	482-022
6	right	482-012	482-132	482-024
1	left	482-001	482-121	482-013
2	left	482-003	482-123	482-015
3	left	482-005	482-125	482-017
4	left	482-007	482-127	482-019
5	left	482-009	482-129	482-021
6	left	482-011	482-131	482-023

1.1.2 Tibial components

Tibia STD, mobile bearing:

Size	Version	Reference number		
		cemented		cementless
		ISO 5832-4 CoCrMo	ISO 5832-4 CoCrMo TiNbN	ISO 5832-4 CoCrMo TPS, BONIT®
1	right	482-098	482-134	482-026
2	right	482-100	482-136	482-028
3	right	482-102	482-138	482-030
4	right	482-104	482-140	482-032
5	right	482-106	482-142	482-034
6	right	482-108	482-144	482-036
1	left	482-097	482-133	482-025
2	left	482-099	482-135	482-027
3	left	482-101	482-137	482-029
4	left	482-103	482-139	482-031
5	left	482-105	482-141	482-033
6	left	482-107	482-143	482-035

Tibia STD, fixed bearing:

Size	Version	Reference number		
		cemented		cementless
		ISO 5832-4 CoCrMo	ISO 5832-4 CoCrMo TiNbN	ISO 5832-4 CoCrMo TPS, BONIT®
1	right	522-902	522-926	522-914
2	right	522-904	522-928	522-916
3	right	522-906	522-930	522-918
4	right	522-908	522-932	522-920
5	right	522-910	522-934	522-922
6	right	522-912	522-936	522-924
1	left	522-901	522-925	522-913
2	left	522-903	522-927	522-915
3	left	522-905	522-929	522-917
4	left	522-907	522-931	522-919
5	left	522-909	522-933	522-921
6	left	522-911	522-935	522-923

Tibia Modular, with locking screw, mobile bearing:

Size	Reference number	
	cemented	
	ISO 5832-4 CoCrMo ISO 5832-3 Ti6Al4V	ISO 5832-4 CoCrMo ISO 5832-3 Ti6Al4V TiNbN
2	482-170	482-202
3	482-171	482-203
4	482-172	482-204
5	482-173	482-205
6	482-174	482-206

Tibia Modular, with locking screw, fixed bearing:

Size	Reference number	
	cemented	
	ISO 5832-4 CoCrMo ISO 5832-3 Ti6Al4V	ISO 5832-4 CoCrMo ISO 5832-3 Ti6Al4V TiNbN
2	522-937	522-942
3	522-938	522-943
4	522-939	522-944
5	522-940	522-945
6	522-941	522-946

Augment:

Size	Height [mm]	Version	Reference number
			cemented
			ISO 5832-3 Ti6Al4V
2	5	left lateral / right medial	482-182
3	5	left lateral / right medial	482-183
4	5	left lateral / right medial	482-184
5	5	left lateral / right medial	482-185
6	5	left lateral / right medial	482-186
2	10	left lateral / right medial	482-187
3	10	left lateral / right medial	482-188
4	10	left lateral / right medial	482-189
5	10	left lateral / right medial	482-190
6	10	left lateral / right medial	482-191
2	5	left medial / right lateral	482-192
3	5	left medial / right lateral	482-193
4	5	left medial / right lateral	482-194
5	5	left medial / right lateral	482-195
6	5	left medial / right lateral	482-196
2	10	left medial / right lateral	482-197
3	10	left medial / right lateral	482-198
4	10	left medial / right lateral	482-199
5	10	left medial / right lateral	482-200
6	10	left medial / right lateral	482-201

Stem extension:

Length	Version	Reference number
		cemented
		ISO 5832-3 Ti6Al4V
40	straight	482-175
75	straight	482-176
110	straight	482-177
110	with offset	482-178

Centraliser:

Size	Reference number
	cemented
	ISO 5834-2 UHMWPE
Ø10	482-179
Ø12	482-180
Ø14	482-181

Insert STD, mobile bearing:

Size	Height [mm]	Reference number		Size	Height [mm]	Reference number	
		ISO 5834-2 UHMWPE ISO 5832-1 1.4441	Vitamin-E XL-PE ISO 5832-1 1.4441			ISO 5834-2 UHMWPE ISO 5832-1 1.4441	Vitamin-E XL-PE ISO 5832-1 1.4441
1	10	482-054	482-301	4	10	482-078	482-310
1	12	482-055	482-302	4	12	482-079	482-311
1	14	482-056	482-303	4	14	482-080	482-312
2	10	482-062	482-304	5	10	482-086	482-313
2	12	482-063	482-305	5	12	482-087	482-314
2	14	482-064	482-306	5	14	482-088	482-315
3	10	482-070	482-307	6	10	482-094	482-316
3	12	482-071	482-308	6	12	482-095	482-317
3	14	482-072	482-309	6	14	482-096	482-318

Insert UC, mobile bearing:

Size	Height [mm]	Reference number		Size	Height [mm]	Reference number	
		ISO 5834-2 UHMWPE ISO 5832-1 1.4441	Vitamin-E XL-PE ISO 5832-1 1.4441			ISO 5834-2 UHMWPE ISO 5832-1 1.4441	Vitamin-E XL-PE ISO 5832-1 1.4441
1	10	482-049	482-319	4	10	482-073	482-334
1	12	482-050	482-320	4	12	482-074	482-335
1	14	482-051	482-321	4	14	482-075	482-336
1	17	482-052	482-322	4	17	482-076	482-337
1	20	482-053	482-323	4	20	482-077	482-338
2	10	482-057	482-324	5	10	482-081	482-339
2	12	482-058	482-325	5	12	482-082	482-340
2	14	482-059	482-326	5	14	482-083	482-341
2	17	482-060	482-327	5	17	482-084	482-342
2	20	482-061	482-328	5	20	482-085	482-343
3	10	482-065	482-329	6	10	482-089	482-344
3	12	482-066	482-330	6	12	482-090	482-345
3	14	482-067	482-331	6	14	482-091	482-346
3	17	482-068	482-332	6	17	482-092	482-347
3	20	482-069	482-333	6	20	482-093	482-348

Insert STD, fixed bearing:

Size	Height [mm]	Reference number		Size	Height [mm]	Reference number	
		ISO 5834-2 UHMWPE	Vitamin-E XL-PE			ISO 5834-2 UHMWPE	Vitamin-E XL-PE
1	10	522-947	522-801	4	10	522-962	522-816
1	12	522-948	522-802	4	12	522-963	522-817
1	14	522-949	522-803	4	14	522-964	522-818
1	17	522-950	522-804	4	17	522-965	522-819
1	20	522-951	522-805	4	20	522-966	522-820
2	10	522-952	522-806	5	10	522-967	522-821
2	12	522-953	522-807	5	12	522-968	522-822
2	14	522-954	522-808	5	14	522-969	522-823
2	17	522-955	522-809	5	17	522-970	522-824
2	20	522-956	522-810	5	20	522-971	522-825
3	10	522-957	522-811	6	10	522-972	522-826
3	12	522-958	522-812	6	12	522-973	522-827
3	14	522-959	522-813	6	14	522-974	522-828
3	17	522-960	522-814	6	17	522-975	522-829
3	20	522-961	522-815	6	20	522-976	522-830

Insert UC, fixed bearing:

Size	Height [mm]	Reference number	Size	Height [mm]	Reference number
		Vitamin-E XL-PE			Vitamin-E XL-PE
1	10	522-831	4	10	522-846
1	12	522-832	4	12	522-847
1	14	522-833	4	14	522-848
1	17	522-834	4	17	522-849
1	20	522-835	4	20	522-850
2	10	522-836	5	10	522-851
2	12	522-837	5	12	522-852
2	14	522-838	5	14	522-853
2	17	522-839	5	17	522-854
2	20	522-840	5	20	522-855
3	10	522-841	6	10	522-856
3	12	522-842	6	12	522-857
3	14	522-843	6	14	522-858
3	17	522-844	6	17	522-859
3	20	522-845	6	20	522-860

1.1.3 Patellar components

Patella STD, mobile bearing:

Size	Reference number		
	cemented		cementless
	ISO 5832-12 CoCrMo ISO 5834-2 UHMWPE ISO 5832-1 1.4441	ISO 5832-12 CoCrMo TiNbN ISO 5834-2 UHMWPE ISO 5832-1 1.4441	ISO 5832-12 CoCrMo TPS, BONIT® ISO 5834-2 UHMWPE ISO 5832-1 1.4441
1	482-109	482-145	482-115
2	482-110	482-146	482-116
3	482-111	482-147	482-117
4	482-112	482-148	482-118
5	482-113	482-149	482-119
6	482-114	482-150	482-120

Patella FIX, fixed bearing:

Size	Reference number
	cemented
	ISO 5834-2 UHMWPE ISO 5832-1 1.4441
1	482-163
2	482-164
3	482-165
4	482-166
5	482-167
6	482-168

1.2 Overview of Instruments

Only the instruments of OHST Medizintechnik AG listed below are to be used for implantation:

Designation	Reference number
ZEN size 1 & 6 instrument set	482-541
ZEN Fixed Bearing size 1 & 6 instrument set	522-120
ZEN Femur instrument set	482-542
ZEN General instrument set	482-543
ZEN Fixed Bearing General instrument set	522-121
ZEN Trial Components instrument set	482-544
ZEN Fixed Bearing Trial Components instrument set	522-122
ZEN Patella instrument set	482-545
ZEN Tibia instrument set	482-546
ZEN Tibia Modular instrument set I	521-103
ZEN Tibia Modular instrument set II	521-104

Designation	Reference number
ZEN Anterior Cut First instrument set	482-742
ZEN Optional instrument set	482-775

1.3 Miscellaneous Accessories

Designation	Reference number
ZEN Femur + Insert X-ray templates	482-507
ZEN Femur X-ray templates	482-559
ZEN Tibia X-ray templates	482-508
ZEN Tibia Fixed Bearing with Insert X-ray templates	522-501
ZEN Tibia Modular X-ray templates	482-738
ZEN Tibia Modular Fixed Bearing with Insert and Stem Extension Straight X-ray templates	522-502
ZEN Tibia Modular Offset Stem X-ray templates	482-739
ZEN Tibia Modular Fixed Bearing with Insert and Stem Extension with offset X-ray templates	522-503
ZEN Knee System Mobile Bearing surgical technique	50000434
ZEN Tibia Modular Mobile Bearing surgical technique	50000435
ZEN Knee System Fixed Bearing surgical technique	50000656
ZEN Tibia Modular Fixed Bearing surgical technique	50000657
Implant passport	50000572

2. Handling

2.1 General information

This implant is part of a system and must only be used with the corresponding original system components.

Only the above-mentioned instruments pertaining to the respective system are to be used for implantation. Before the instruments are used, the corresponding instructions for use (50000354) must be followed.

Caution: implants must always be stored in their full, unopened protective packaging. The implant packaging must not be exposed to direct sunlight. Prior to insertion of the implant, the packaging should be inspected for damage, as this could affect sterility.

When removing the implant from the packaging, the user must check that the implant corresponds to the designation on the packaging (Art. no. / serial number / size).

The respective size and, in the case of inserts, the respective height are shown by a colour code on the label attached to the outside of the implant's cardboard packaging. Further information on colour coding can be found in the respective surgical technique of the knee system (see section 1.3. Miscellaneous Accessories). Compliance with the relevant hygiene regulations is required during removal of the implant from the packaging. Care must be taken to protect all implant surfaces against damage, since this could be decisive for possible failure. Therefore the prosthesis must not come into contact with any objects which could damage its surface. Prior to use, every implant must be visually inspected for damage.

Modifying an implant can not only reduce its service life, but can also lead to immediate or subsequent failure of the prosthesis under stress. The implant must therefore not be modified either mechanically or otherwise.

Implants from damaged packaging, unsterile, contaminated, damaged or carelessly handled implants or implants subjected to unauthorised modification must not be used.

Caution: the implants are intended for single use only! The individual loads on the functional surfaces of a patient affect these surfaces to such an extent that re-use of the implants is excluded. The marks of stress on the functional surfaces cannot be reliably detected by visual inspection alone. After an explantation, it therefore has to be assumed that there is damage precluding re-use.

For implant components that are designated for only one side of the body, the corresponding side is marked on the implants with "L" for the left knee and "R" for the right knee. The orientation of the implants must always correspond to the side of the body of the knee to be replaced. Implant components on which no particular side of the body is marked can be used for either the left or the right knee.

2.2 Authorised Component Combinations

Together with its instruments, the ZEN Knee System is designed as a complete system and does not allow for the substitution of components from other systems of manufacturers. Only the above-mentioned system components may be used for implantation.

Compatibility instructions ZEN Mobile Bearing:

Each of the femoral components must be combined with inserts of the same size.

The inserts can be combined with tibial components of the same and of the next smaller size, as well as with all larger tibial components. This applies to both the STD inserts and the UC inserts.

Compatibility instructions ZEN Fixed Bearing:

Each of the tibial components must be combined with inserts of the same size.

The inserts can be combined with femoral components of the same size as well as with the next larger and the next smaller size. This applies to both the STD inserts and the UC inserts.

General compatibility instructions:

The femoral components, stem extensions and augments can be used for both the Fixed Bearing variant and for the Mobile Bearing variant. Each of the patellar components are to be used with the femoral components of the same size.

We guarantee compatibility of our products only in combination with our own CE marked products and with the products we have approved for combined use and which have been authorised accordingly. In this regard, please note the instructions for use provided by the endoprosthesis manufacturers and the combination matrix approved by OHST.

Due to reasons relating to product safety and product liability, it is prohibited to use implants manufactured by OHST Medizintechnik AG in combination with components of other manufacturers that have not been approved by OHST.

2.3 Instructions for Use

Before the instruments are used, the enclosed product information and instructions for use must be observed.

The surfaces with a porous coating (TPS, BONIT®) and the roughened implant surfaces must not come in contact with clothing or other fibre-shedding materials.

Caution: whenever possible, avoid contact with the BONIT®-coated implant parts. These areas should be touched only with powder-free latex gloves.

Follow the manufacturer's instructions on using the bone cement to ensure correct performance of the cementing step. The use of bone cement mixed under vacuum is recommended.

Prior to introducing the cement (for cemented fixation) or inserting the implant (for cementless fixation) the implant bed must be adequately rinsed. Care must be taken that all loose particles (e. g. bone splinters, abrasion particles of the tools etc.) are removed from the prepared implant bed. In the case of cemented fixation, the implants must be inserted into the cement bed centrally and upright. Following the cementing step, all protruding or loose cement particles must be removed from the wound area.

Caution: when using high-frequency surgical instruments (e.g. cautery knife), it is essential to ensure that they do not come into contact with the implants or instruments. This can cause such severe damage to the implants or instruments that failure (e.g. fracture) may ensue. If an implant has been damaged, it must not remain in the patient but must be replaced by a new, intact implant. Damaged instruments may only continue to be used if they can still perform their intended purpose without compromise.

2.4 Surgical Technique

2.4.1 ZEN Knee

Approach

- Parapatellar incision with medial arthrotomy or another approach, e.g. for a revision procedure or marked valgus malposition

Flexion gap - tibia processing step

When using the ZEN Tibia Modular with a stem extension and optional augments, the information provided in section 2.4.2 ZEN Tibia Modular applies.

- Mounting and positioning of the tibial goniometer
- Adjustment of rotation in the frontal plane
- Adjustment of the slope
- Adjustment of the cutting height
- Attachment of the tibial cutting block
- Performance of the tibial cut
- Verification of the tibial cut

Flexion gap - femur processing step

- Determination of the femur size
- Intramedullary opening
- Align and fix femoral cutting block in a.-p. position
- Anterior and posterior femoral resection

Caution: This step is a temporary step which has to be finalised in a later work step.

- Review of the flexion gap

Extension gap - distal femur processing step

- Distal positioning of the cutting block
- Review of the angle and distal cut
- Review of the extension gap

Finalisation of the femur

- Attachment of the femoral cutting block finisher
- Execution of the cuts

Caution: It is important to carry out all cuts precisely, as the previous A/P incision was only temporary.

- Verification of the cuts
- Preparation of the trochlear groove

Tibia processing step

When using the ZEN Tibia Modular with a stem extension and optional augments, the information provided in section 2.4.2 ZEN Tibia Modular applies.

- Selection of the plateau
- Tibia finalisation

Patella processing step

For reduced bone thickness, the thinner full-poly patella should be employed, as it guarantees greater preservation of the bone. Prosthetic patellar replacement is contraindicated where the patella is too thin or in the case of a flat patella.

- Removal of the osteophytes
- Align and fix patellar resection guide perpendicular to the patellar tendon
- Measurement of the patellar thickness (at least 13 mm of bone must remain after deducting the thickness of the patellar component)

Components	Height of patella STD	Height of patella FIX
Size 1	11.2 mm	9.7 mm
Size 2	11.2 mm	9.7 mm
Size 3	11.2 mm	9.7 mm
Size 4	11.2 mm	9.7 mm
Size 5	11.7 mm	10.2 mm
Size 6	11.7 mm	10.2 mm

- Adjustment of the height of the drilling gauge and insertion of the pins
- Removal of the patellar resection guide
- Positioning of the patellar cutting block
- Performance of the cut
- Alignment of the base plate of the trial patella

- Clamping of the patellar drilling guide into the patella forceps
- Positioning of the patella base plate and compression of the patellar forceps
- Drilling of the pin holes
- Removal of the patellar forceps and drilling guide
- Insertion of the trial patella insert

Attachment of the trial components

When using the ZEN Tibia Modular with a stem extension and optional augments, the information provided in section 2.4.2 ZEN Tibia Modular applies.

- Components are inserted in the following order: Trial tibia, trial insert, trial femur and, if applicable, the trial patella

Caution: The implant components must not be used together in any random combination. Compatibility instructions are provided in section 2.2.

Caution: Should the height of the insert have to be changed or the trial components removed after successfully checking the overall mobility of the joint, the components are removed in the reverse order of the insertion. In other words: first remove the trial femur, followed by the trial insert.

Implantation

When using the ZEN Tibia Modular with a stem extension and optional augments, the information provided in section 2.4.2 ZEN Tibia Modular applies.

- Components are inserted in the following order: Tibial components, insert, femoral component and, if applicable, the patellar component

Caution: The UC (ultra-congruent) inserts are to be used only after the resection of both cruciate ligaments. They do not provide any room for a posterior cruciate ligament that is still intact.

2.4.2 ZEN Tibia Modular

Access

- Parapatellar incision with medial arthrotomy or another form of access, e.g. for a revision procedure or marked valgus malposition

Flexion gap - tibia processing step

Intramedullary alignment:

- Positioning of the intramedullary tibial alignment guide
- Medial-lateral centred opening of the tibia, aligned 1/3 anterior and 2/3 posterior
- Use of the medullary drills in ascending order, starting with the smallest diameter
- Mounting of the intramedullary tibial alignment gauge on the in-situ medullary drill
- Adjusting the rotation in the frontal plane
- Performing varus-valgus adjustment
- Fixing the rotation and varus-valgus adjustment
- Adjustment of the cutting height

- Checking the incision
- Attaching the tibia modular cutting block
- Removing the IM alignment guide and the medullary drill
- Checking the alignment
- Performing the tibial incision
- Checking the tibial incision

Tibia processing step:

- Selection of the plateau
- Positioning and fixing of the trial tibia modular without/with offset
- Drilling the tibial shaft
- Drilling the shaft extension without/with offset
- Driving in the trial cone

Attachment of the trial components

- Sequence: Trial insert, trial femur (the trial tibia and possibly the trial patella had already been inserted in the previous step)
- Selection of the trial insert and trial femur according to the previously determined femur size.

Caution: The implant components must not be used together in a random combination. Information on compatibility can be found in section 2.2.

- Combine trial inserts for size and height and insert them into trial tibia
- Place the trial femur with the femoral impactor on the femur, which is bent by at least 90°
- Axis check
- Insert trial patella if necessary
- Check the overall mobility

Caution: Should the height of the insert have to be changed or the trial components removed after the overall mobility of the joint has been successfully checked, then remove the components in the reverse order to their insertion. In other words, first remove the trial femur and then the trial insert.

Implantation

- The components are inserted in the following order: Tibial component with shaft extension and optional augments and centralisers, the insert, the femoral component and - where applicable - the patellar component
- The shaft extension is placed onto the tibia and fixed in place with the screw provided and then tightened with the torque screwdriver
- The alignment of the shaft extension with offset on the tibia is performed using the mounting device
- With shaft extension L110, the corresponding centraliser is also mounted
- Optional augments are attached using the torque screwdriver

Caution: The UC (ultra-congruent) inserts are used only following the resection of both cruciate ligaments. They do not provide any space for a posterior cruciate ligament that is still intact.

2.4.3 General information

It is important to ensure that the surfaces are not damaged. Appropriate impactors are available for impacting the implants. Position the impactor at the centre and align it parallel to the mechanical axis. The patella is pushed into the bone using the patellar forceps and the pressing aid.

In order to prevent excessive heating of the bone, all resections must be performed while cooling with physiological saline solution. The thickness of the saw blades must be adjusted to the instruments used!

The smaller height of the fixed-bearing patellar component may provide greater bone preservation.

For further information see the surgical technique of the knee system (see 1.3. Miscellaneous accessories).

3. Packaging and Sterility

Depending on the sterilisation method used, implants are packaged in a triple transparent bag made of plastic laminated film (sterilisation by irradiation at least 25 kGy) or in a double transparent bag made of Tyvek® (ethylene oxide sterilisation) inside a carton. The instruments are supplied in a non-sterile condition in protective packaging. They must be cleaned and sterilised prior to use in accordance with the respective instructions for use (50000354). The stated use-by date presumes that the packaging is intact and unopened and that the product was stored under appropriate conditions.

Caution: the implants must not be re-sterilised. Reprocessing of components which have not been implanted and whose packaging has been opened may only be carried out by the manufacturer, as individual validated processes must be repeated.

The outer bag of the triple transparent bag packaging is to be removed by the non-sterile personnel together with the carton. For the double transparent bag packaging, only the carton is to be removed by the non-sterile personnel. The second bag must be opened in such a way that the sterility of the inner bag is not compromised. The inner bag is removed and opened by the sterile personnel. The implant must then be presented to the surgeon, who can then directly remove the sterile implant.

4. Preoperative Planning and Postoperative care

The preoperative planning based on X-rays, CT data and other information is indispensable and provides essential information concerning the suitable prosthesis, its positioning and the possible component combinations, thereby facilitating preselection of the implant size.

Surgery may only be performed once it has been established that the patient is able to tolerate the implant material.

Both conventional and digital X-ray templates (see 1.3. Miscellaneous Accessories) are available for planning the surgery. To check the correct fit, trial prostheses and additional implant components should be kept handy in the event that other sizes are required, or where it is not possible to use the planned implant.

Recognised procedures must be used for postoperative care.

5. Indications

Total prosthesis is indicated in patients who have lost their mobility for any of the following reasons:

- Degenerative, post-traumatic or rheumatoid arthritis
- Avascular necrosis of the femoral condyle
- Post-traumatic loss of the joint configuration, particularly in the case of patellofemoral erosion, dysfunction or prior patellectomy
- Moderate varus, valgus or flexion malpositions

This system may also be indicated for the correction of earlier surgical attempts which have been unsuccessful.

Total knee replacement is not intended to withstand the same levels of activity and loads as normal, healthy bones would; however, in many cases, it can restore mobility whilst simultaneously relieving pain. Total knee replacement should be the course of action only where all other treatment approaches which could preserve the joint, whether surgical or conservative, and which are deemed to be timely and appropriate in medical terms, have failed.

6. Contraindications

- Acute or chronic, local or systemic infection
- Severe muscle, nerve or vascular diseases that endanger the extremity concerned
- Missing bone substance or poor bone quality that threatens the stable fit of the prosthesis
- Hypersensitivity to the materials used
- Loss of the ligamentous apparatus
- Overweight
- Alcoholism or other addictions
- Intense physical activity (e.g. competitive sports or heavy physical labour)
- Any underlying condition that might compromise the function of the implant

7. Risks and conditions that may affect the success of the surgery

General risk factors and conditions:

- Overweight
- Alcoholism or substance abuse
- Patient groups with mental disorders or addictions
- Pregnancy
- Intake of high-dose cortisone or cytotoxicics
- History or risk of infectious diseases with potential joint manifestation
- History of deep vein thrombosis and/or pulmonary embolism
- All general surgical risks

Risk factors and conditions specific to knee arthroplasty:

- Disorders of the bone metabolism (osteoporosis, osteomalacia)
- Circulatory disorders of the affected limb
- Neurological disorders of the affected limb
- Muscular dysfunction of the affected joint
- Growth in children and adolescents
- Expected extreme stresses, e.g. due to work or sports
- Epilepsy or other causes of repeated accidents with an increased risk of fracture

- Joint deformities that complicate fixation of the implant
- Weakening of the supporting structures due to a tumour

Caution: If any of the conditions stated above apply to the patient, the latter must be informed of the likely impact of these factors on the success of the surgery and the service life of the implant.

8. Undesirable Effects

The undesirable effects listed below are among the most typical and commonly occurring consequences of a surgical procedure:

- Infection
- Venous thrombosis and pulmonary embolism
- Cardiovascular disorders
- Haematoma
- Paraesthesia
- Numbness
- Swelling
- Nerve damage
- Oedema

The undesirable effects listed below are amongst the most typical and commonly occurring consequences of total knee replacement surgery:

- Tibial or femoral fracture
- Sub-dislocation or dislocation of the patella
- Instability, migration or loosening of components
- Detachment of components
- Insufficient freedom of movement
- Metal hypersensitivity reactions
- Patellofemoral complications
- Neurovascular complications
- Paralysis/palsy of the peroneal nerve
- Fat embolism
- Arterial insufficiency and damage
- Arthrofibrosis / accretions
- Implant fractures
- Implant noises
- Reduced quality of life (pain, sleep disorders and a restricted range of motion, especially when lying down)
- Inflammations
- Metallosis
- Increase in the level of metal ions in the blood
- Pseudotumours

Caution: Revision surgery may be required if specific undesirable effects occur.

9. Patient Information and Documentation

The serial numbers of the implants used must be documented in the patient's records. The corresponding labels are enclosed in the packaging of the sterile implants for this purpose.

The patient must be informed of the benefits and risks of the procedure. If the implant is deemed to be the best solution for the patient, even though some of the contraindications described above may in part apply to the patient, it is extremely important that the patient is made aware of the effects of these circumstances on the outcome of the surgical procedure and the associated risks. Patients receiving a knee replacement must be informed that the service life of the implant will depend on the patient's weight and how active the patient is in everyday life. The patient must be informed of which activities he/she can undertake in order to reduce the effects of these aggravating circumstances.

All information provided to the patient must be documented in writing by the surgeon performing the surgery. Following surgery, the patient must be given an implant passport containing all the necessary information about the implant. Adhesive labels to document the implant used are enclosed. Undesirable effects that can harm the patient may occur in MRI examinations. Possible effects include artefacts, heating up of the implant, induction of electrical currents, loosening of the implant. Before use, the information provided by the equipment manufacturer should be studied. In case of doubt and as part of an individual risk assessment, the suitability of comparable implants in the respective MRI machine should be checked. The patient is to be informed of the risks.

10. Explanation of the Label Symbols

An explanation of the symbols used by OHST Medizintechnik AG can be found in the annex (p. 143).

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